



# User Guide



# Jogger



**SPECIAL  
TOMATO**

# Jogger™

## Affordable. Maneuverable. Transport

If you're dashing around town, or just going in and out of shops, you really want a lightweight, maneuverable stroller that does not compromise positioning. With it's super-light aluminum frame, the Special Tomato® Jogger folds with ease and can be placed neatly into your vehicle. The comfortable padded seat reclines and there is a protective canopy with a window that allows you to keep an eye on your child. For children who need minimal to moderate positioning, the Soft-Touch® Liner easily attaches. And the Special Tomato® Jogger still folds with the Soft-Touch® Liner in place!

### JOGGER HIGHLIGHTS

- Compact Folding– easy one step folding
- Maneuverable– front wheel swivels
- Comfortable– thick upholstery
- All-Terrain– metro to park trails
- Contoured body support– with Soft-Touch® Liners
- Machine washable fabric– easy to remove

### STANDARD FEATURES

- 3 wheels (12" pneumatic tires)
- Front disc brake
- Foot well with side-guard
- Functional rear suspension
- Reversible Fleece/ Canvas Insert
- Rear parking brake
- Recline- seat back
- Canopy- folds with jogger
- 5-point harness
- Height adjustable push handle
- Thick padded seat & back surface
- Quick release rear wheels
- Folding Latch
- Latex free

### OPTION FEATURES

- Padded Head Support
- Medical Necessities Bag
- Soft Touch® Liner back pad (available in 2 sizes)
- Soft Touch® Liner seat pad (available in 2 sizes)

5632-70	TOMATO JOGGER
5632-71	MEDICAL BAG
5632-72	HEAD SUPPORT

5618-05	5618-11	SIZE 1 SEAT CUSHION
5618-06	5618-12	SIZE 1 BACK CUSHION
5618-07	5618-13	SIZE 2 SEAT CUSHION
5618-08	5618-14	SIZE 2 BACK CUSHION



### SPECIFICATIONS


Weight.....	27lbs
Seat to Foot-well .....	15"
Seat Width .....	14"
Seat Depth .....	12"
Back Height .....	22"
Total Width .....	26"
Max user weight .....	110lbs or 50kg
Total Length .....	42"
Total Height .....	43"
Handle Height to Floor.....	35" – 43"
Folded Dimensions .....	33" x 21" x 16"
(without rear wheels)	
Recline Range .....	100° – 140°



• Height adjustable push handle



• Compact folding



• Functional rear suspension




## Contents

---

2	Frame
3	Front Wheel Attachment
3	Front Wheel In-line Locking
4	Rear Wheel Attachment
4,5	Safety Harness
6	Parking Break
6	Recline Seat
7	Removing Seat Fabric
7,8	Shopping Basket
8	Canopy Attachment
8	Handle Height Adjustment
9	Folding & Optional Liners

**In the interest of your priceless cargo, please read these warnings before using your Tomato Jogger.**


- Keep this booklet for future use.
- Accessories which are not approved by the manufacturer or distributor should not be used.

-  • **WARNING:** It may be dangerous to leave child unattended.
-  • **WARNING:** Attaching objects to the handle affects the stability of the jogger.
-  • **WARNING:** To Prevent serious injury, always use the 5-point safety harness.
  - Always use the wrist strap.
  - Child should remain seated in the jogger, do not use as a step stool.
  - The Special Tomato™ Jogger seats one passenger.
  - Never use Roller Blades while using the Special Tomato™ Jogger.
  - Total weight load recommended for this Jogger is 110 lbs or 50 kg.

## Frame

---

Unpack Jogger from box and set wheels aside.



1. Push foot plate of Jogger away from handle by pushing curved crossbar up and out.






2. Audible click will be heard when Jogger frame is completely unfolded.



3. Flip Jogger to upright position.



4. Clip lower belting around frame to secure seat fabric.




## Front Wheel

To install front wheel, lay jogger on it's back so handle and rear axles are touching the floor.

5a. Insert the Front Wheel into the fork tips by aligning the axle flats with the fork tip flat surfaces.

5b. Tighten by turning the wheel lever and securing the other side.

5c. Close the wheel lever making sure that it is flush with the tubing. Be sure not to over tighten.



Remove Black Caps

## Front Wheel In-line Locking

7. Behind the front wheel pivot point is where the in-line locking mechanism is located. **Pull and rotate knob to lock and unlock.** Image shows unlocked.

8. Image now shows the wheel locked in-line. Pull and rotate mechanism to unlock.

## Rear Wheels



9. To install, gently slide the rear wheel's axle into the frames supporting axle hole until you hear it click onto place.




10. To remove, simply push down the spring lock tab while pulling out on the wheel.



## Safety Harness

Adjusting safety harness:



11a. Unclip hooks.

11b. Remove fleece lining (Velcro tabs shown above).



12

**12. Unthread harness straps from fleece lining.**



13

**13a. Adjust shoulder straps to desired height using the loops on the seat back.**

**13b. Slide the plastic guide up and down to adjust strap length.**

**Once fitting is complete reassemble the fleece lining .**



14


**14a. Clip both shoulder straps to the eye holes located on the buckle mechanism.**

**14b. To adjust the length of the center strap, move its plastic guide up and down.**

## Rear Wheel Parking Brake


To apply the parking brake:

15. Press down on the foot bar.



Disengage the brake:

16. Pull up on the foot bar.



## Reclining Seat


17. Simply adjust the recline strap length to desired degree of recline.



18. Unlock the buckle shown in Figure 17 for maximum recline.



## Removing Seat Fabric




19. Using a phillips head screw driver unscrew hardware from the frame.



20. Unsnap fasteners on both sides of jogger.


20a. Unlock seat strap from frame, buckle is suspended across the back of the jogger.



21. Under the Foot Well , Unloop Velcro Strap through the D-ring.

**NOTE: REMOVE ALL RIGID PARTS BEFORE WASHING**  
Joggers are shipped with Seat and Basket pre-assembled. Assembly/ Removal instructions are given if any pre-assembly items are removed for cleaning or replacement. "SEE PAGE 10"


## Shopping Basket



22. Secure front basket flap by snapping around the lower foot well crossbar.

(Detail of front shopping pocket snaps)






23. Wrap rear basket straps around the frame and loop back through the D-ring. Secure the Velcro strap.

*Note: The D-ring loops are at the rear of the basket.*




## Canopy

24. Fasten the canopy to the portion of the upper frame that is not covered by fabric. Small holes in the crest of the frame tubing accept the canopy. Be sure there are NO obstructions (fabric or brake lines) between the canopy clip and the frame. There is a small plastic nipple on the canopy mechanism that fits into the hole in the frame.



25. Hook the wire on the back side of the canopy mechanism and then close the lever by pressing it to the frame (Figure 25). **DO NOT FORCE TO CLOSE!** Double check for obstructions (fabric or break line/ be sure nipple is inserted to frame).

## Handle Height



26. Unlock lever on both sides of handle. Adjust to height. Then relock.





## Folding

27. Squeeze mechanism located below the handle. (Mechanism located on Both Sides)



Fold jogger back. Red clip secures the closure.



## Optional Special Tomato Liners

28. Back Liner:

Using the provided straps that come with the back liner – Feed straps through the slots in the jogger seat back. [Red Arrows] Extra slots provide adjustment. (Do Not feed Jogger shoulder straps through liner.)

29. Seat Liner:

Use attachment straps provided with the seat liner to fasten the optional back pad to the jogger. [Red Arrows]

Then feed the original jogger seat buckles up through the seat liner. [White Buckles]

## Wrist Strap

- Make sure safety wrist strap is secured properly around your wrist whenever your baby is in the jogger.

## Washing

- You can machine wash your jogger's seat in cold water with a mild detergent. To minimize shrinkage, simply reattach the seat to the frame to dry. ( DO NOT put seat in the dryer)
- Be sure to remove the rigid back board from the seat before washing. Unzip seat back and remove rigid back support.
- You may clean the jogger frame using a soft, dry cloth.

## Wheel Care

- Moisture can rust the wheel bearing case. Spray regularly with a silicone lubricant to keep your wheels looking clean and functioning like new.
- Also, if you trek a massive amount of miles, your jogger will benefit from a visit to the local bike shop for a tune up. Refer to the number on the tire wall for optimal tire pressure.

## Storing

- Sunshine can fade seat fabric and dry out tires, so storing your jogger indoors will prolong color vibrance.



**WARNING:** Prevent serious injury or death from falls or sliding out

- Always secure child in the restraint
- Never leave child unattended
- Always keep child in view
- Max User Weight 110lbs or 50kg



**Bergeron**  
by design


Bergeron By Design  
15 South Second St.  
Dolgeville, NY 13329  
1.866.529.8407

# Thank You!

We wish you happy trails with your new Special Tomato™ Jogger! Your jogger is a great product, and we are proud to have it in our line. If you have any questions or comments, please give us a call toll free at.

1.866.529.8407

---



Special Tomato is a Product Line from Bergeron By Design

15 South Second Street • Dolgeville, NY 13329  
T: 1.866.529.8407 • F: 1.315.429.8862

