

Invacare® | **Leo™**



Safe, solid, secure: the *Invacare Leo*

The *Invacare Leo* is a 4-wheel scooter designed for all those that value their independence and wish to get out and about unaided. Safety is a key feature of the **Leo**, but this does not detract from its stylish and sporty looks. **Leo** offers people the freedom and confidence to enjoy their essential daily outings and leisure excursions. With high standards of design and engineering safety, the **Leo** undoubtedly sets a new standard.



Yes, you can.®



For a safe, powerful drive

The *Invacare Leo* offers a smooth, safe drive in a scooter that is easy to control and manoeuvre. The **Leo** intelligent lighting package includes indicators that automatically switch off and a rear brake light for safe driving. The lighting system is designed so that if a bulb should blow, all other lights still continue to function. In addition to a visual warning for the battery level display, drivers are also alerted with an audible warning (beep) when the battery needs recharging.

Compact

The **Leo** is a compact mini-scooter that can be easily dismantled for transportation during longer trips or weekend breaks. With a smooth and manoeuvrable drive, getting round tight corners in the local shopping centre, museum or leisure complex becomes refreshingly easy. Access to the battery and electronics for maintenance or servicing is much easier with the **Leo** as only one part needs to be removed.

Reliable

Customers can trust the **Leo** to take them wherever they need to go. From a trip to the shops or a ride to the park, the **Leo** is rendered sturdy and reliable thanks to its exceptionally high quality electronic engineering and power. These elements are protected from corrosion and water,

(splash protection) which safeguard the vital workings of the scooter against dirt and water. A two-step disengaging lever protects the scooter from free-wheeling if inadvertently knocked out of drive mode.

Features and Accessories



Finite tiller movement

Customers can easily adjust the tiller to suit their needs.



Brake light

Rear brake light warns other people that the scooter is slowing down, thus preventing collisions.



Comfort seating

Swivel function and sliding rails ensure an optimal seating position.



Two-step disengaging lever

This unique concept prevents the scooter from accidentally freewheeling.

Features and Accessories



Built in anti-splash guards
Protects electronics and the transaxle from dirt and water.



Disassembly into five parts
No tools are required to dismantle
Invacare Leo



Crutch & cane holder
Keeps your crutch or cane stable during the drive.



Lockable front box
Secure box for personal valuables.



Seatbelt
Keeps a person well-positioned in the seat.



Large mirror
Can be mounted left or right.



Heavy-duty storage cover
Keeps the scooter clean and dry.









Charging off-board
Simply remove batteries to charge with the off-board charger kit.



Seat suspension module
Adds extra comfort during driving.

Technical data

						
Invacare Leo	470 mm	410 mm	460 - 560 mm	475 mm	590 mm	1220 mm


					
Invacare Leo 59 kg (excl. batteries) 83 kg (incl. batteries)	136 kg	2 x 36 Ah	2620 mm	60 mm	70 mm

				
Invacare Leo	10°	240 W	36 km	8 km/h

Colour of shrouds



Upholstery



INVACARE Limited
Pencoed Technology Park - Pencoed -
Bridgend CF35 5AQ
Tel: +44 1 656 776 222
Fax: +44 1656 776 220
E-mail: uk@invacare.com
www.invacare.co.uk
Sales Order E-mail:
ordersuk@invacare.com

INVACARE Ireland Ltd
Unit 5 - Seatown Business Campus
Seatown Road - SWORDS - County
Dublin
Tel. +353 1 8107084
Fax +353 1 8107085
Email: Ireland@invacare.com
www.invacare.ie

© 2011 Invacare International Sàrl
All rights reserved.

All information quoted is believed to be
correct at time of print.
Invacare reserves the right to alter
product specifications without prior
consultation.

Leo - UK - 07/2011 - 1503875

