

Invacare® | ***Birdie® and Birdie® Compact***



5 year warranty*



Transportability and patient comfort go hand-in-hand

With maximum space for the user, both the *Invacare Birdie* and the *Birdie Compact* offer a truly comfortable lift and transfer to or from beds, chairs or even the floor. The lifters are designed to ensure that folding and unfolding can be carried out easily and without the need for tools. In addition, the lifter can be dismantled into two parts if needed – also without tools.



Yes, you can.®



Total user comfort

The high boom maximizes space for the client even when the hoist is in the highest position. This gives an ergonomic lift and transfer as the client can be easily rotated in the hoist. This makes placing the client into a bed/wheelchair much easier for the carer. The optimization of space creates a reduction in risk of injury as the client's knees will not be close to the actuator of the hoist.

Ease of handling: simple to fold, unfold and dismantle without tools

Handling the lifter has been made as smooth and easy as possible as the hoist can easily be folded and unfolded without the need for tools. The **Birdie** takes up minimal space when stored and is easy to push and transport. If required, the mast and boom can be easily detached from the base (no tools required).

TÜV-certification for maximum security

For total security, the hoists come with a snap-lock system, thus minimizing the risk of accident. The lifters are CE-certified in line with the Medical Device Directive 93/42/ECC. To ensure

the hoists comply with EN ISO 10535 (the European standard for hoists) the hoist has been rigorously tested both by Invacare and the German TÜV institute, an accredited external body.

Features and Options




Long boom – high user comfort

The user can easily be rotated and positioned. No risk of injuries to the knees.



Easy to fold and dismantle

Ease of transport and storage. No tools required when folding / unfolding or when dismantling the hoist.



Snap lock

Optimal security to avoid risk of accidents.



Emergency stop and emergency lowering

Manual or electric emergency lowering (compact version features only electric lowering). Available with fixed or detachable battery. Off board battery charger is available.

Features and Options



Click on spreader bar
Easy to change the spreader bar.



Manual leg spread
Simply step on the pedal to adjust the leg spread.



Handle for leg spread
Handle can easily be attached to operate the manual leg spread.



Electrical leg spread
Available for Invacare **Birdie** (option).



2 point quick release spreader bars
Easy to attach and detach. Available widths: 350 mm, 450 mm, 550 mm.



Push bar and hand control
Superior ergonomics for carers.



Castors
100 mm or 75 mm castors available. (braking castors are shown)













Scale
If required, a scale is available approved. According to the medical standard EN 45501, certified and sealed in order to be used to medicate.



4 point quick release spreader bars
Easy to attach and detach. Available widths: 450 mm and 550 mm.

Technical data

					
Birdie	1150 mm	655 mm	42 kg	465 - 1685 mm	100 mm
Birdie Compact	1000 mm	530 mm	36 kg	515 - 1585 mm	100 mm

					
Birdie	1040 mm	170 kg	40 cycles	2,9 Ah	1400 mm
Birdie Compact	890 mm	150 kg	40 cycles	2,9 Ah	1070 mm

All measures are taken with 75 mm castors and 2 point spreaderbar.
For 100 mm castors add 15 mm for height measures and 20 mm for width.

Frame colour



Light grey /
Dark grey

*Warranty:

Actuator: 5 years
Frame: 3 years
Battery: 1 year

INVACARE® Limited
Pencoed Technology Park - Pencoed -
Bridgend CF35 5HZ
Tel: +44 1 656 776 222
Fax: +44 1656 776 220
E-mail: uk@invacare.com
www.invacare.co.uk
Sales Order E-mail:
ordersuk@invacare.com

INVACARE® Ireland Ltd
Unit 5 - Seatown Business Campus
Seatown Road - SWORDS - County
Dublin
Tel. +353 1 8107084
Fax +353 1 8107085
Email: Ireland@invacare.com
www.invacare.ie

© 2011 Invacare® International Sarl
All rights reserved.

All information quoted is believed to be
correct at time of print.
Invacare® reserves the right to alter
product specifications without prior
consultation.

Birdie - UK - 03/2011 - 1504431