

Invacare® | **Pronto® M41™**



New versatile mobility

The Invacare **Pronto M41** is the latest addition to the **Pronto** range of power wheelchairs. Combining mid-wheel drive and **SureStep®** suspension allows the **Pronto M41** to be a compact, manoeuvrable indoor chair with excellent outdoor driving capabilities.

Its sleek, stylish, candy red shrouding ensures the **Pronto M41** looks as engaging as its performance. The cross compatibility of a number of options and accessories means the **Pronto M41** fits straight into the expansive Invacare range of mobility products.



Yes, you can.®



Optimum Comfort

Pronto M41 has been designed with SureStep technology to ensure excellent stability and manoeuvrability in the narrowest of places. The tension adjustable backrest is a standard feature and is angle adjustable -10° to 30°. The ability to swivel 90° in either direction allows easy access and this combined

with height adjustable armrests, a one piece flip-up footboard and a standard 50 mm cushion ensures maximum comfort for day long use. Suitable for indoor and outdoor environments, the **Pronto M41** provides a perfect solution for those clients who are looking for a chair that is extremely comfortable and versatile.

Features and Options



Center wheel drive

Highly manoeuvrable with centre wheel drive.



SureStep®

Maximum stability and safety with the **SureStep** feature.



Tension adjustable backrest

Velcro adjustable backrest for individual support.

Technical data

| | | | | | | |
|-------------------|---------------|---------------------------|---------------------------|---------------|--------------------------|---------------------------|
| | | | | | | |
| Pronto M41 | 440 - 480 mm | 410, 460, 510 mm | -10° - 30° ^{4,5} | 940 - 1100 mm | 480, 540 mm ¹ | 275 - 465 mm ² |
| | | | | | | |
| Pronto M41 | 610 mm | 980 - 1030 front riggings | 1020 mm | 50 mm | 7° / 12.3% ³ | |
| | | | | | | |
| Pronto M41 | Fixed tilt 5° | 2 x 32 Ah | 2 x 250 W | 21 km | 6 km/h | |
| | | | | | | |
| Pronto M41 | 82 kg | 130kg | | | | |

Frame colour



Candy red

Invacare Limited
Pencoed Technology Park - Pencoed -
Bridgend CF35 5AQ - United Kingdom
Tel: +44 1 656 776 222
Fax: +44 1656 776 220
E-mail: uk@invacare.com
Sales Order E-mail:
ordersuk@invacare.com
www.invacare.co.uk

Invacare Ireland Ltd
Unit 5 - Seatown Business Campus
Seatown Road - SWORDS - County
Dublin - Ireland
Tel. +353 1 8107084
Fax +353 1 8107085
Email: Ireland@invacare.com
www.invacare.ie

© 2013 Invacare International Sàrl
All rights reserved.

All information quoted is believed to be correct at time of print.
Invacare reserves the right to alter product specifications without prior consultation.

Pronto M41 - UK - 03/2013 - 1552649



- 1/ measured without seat cushion
2/ certain legrest lengths are restricted to a special legrest type
3/ the max. safe slope depends on the wheelchair configuration.
For further details please see the user manual
4/ manual settable
5/ Fixed tilt position 5°

