

Invacare® | **Phoenix® Plus**



A great deal more for less

Invacare has taken the original *Invacare Phoenix* and re-styled and updated it to produce a more stylish, more flexible and comfortable power chair, while retaining the core features of easy to use controls, adjustability and quick and easy dismantling. Last but not least, the *Invacare Phoenix Plus* is now available at an even more competitive price.



Yes, you can.®



The smarter choice

A practical, compact power wheelchair, the **Invacare Phoenix Plus** is suitable for both indoors and outdoors. The **Phoenix Plus** comes with 40 Ah batteries as standard, increasing your mobility and independence. The additional features such as angle adjustable footplates, padded clothes guards and armrests improve all-round comfort and positioning, minimising fatigue and tiredness.

Practicality

The **Phoenix Plus** is the ideal choice for manoeuvrability and versatility with its compact footprint.

This practical wheelchair can be easily dismantled for transportation or storage.

Features and Options



Folding Back

Reduces overall size when folded for transport or storage.



Height Adjustable Armrest

The armrests are easy to remove and adjustable in height.



Alloy wheels

Durable alloy rims to ensure a trouble-free ride.

Technical data

| | | | | | | | |
|--|--------------|----------|--------------------|--------------------|-----------------|-------------|---------------------|
| | | | | | | | |
| Phoenix Plus 400 / 440 / 480 mm | 400 - 440 mm | 500 mm | 400 - 460 mm | 220 - 320 mm | 360 - 460 mm | 12° (fixed) | 570 mm ¹ |
| | | | | | | | |
| Phoenix Plus 970 / 1230 mm | 1060 mm | 730 mm | 68 kg ² | 114 kg | 2 x 12V / 40 Ah | | |
| | | | | | | | |
| Phoenix Plus 1940 mm | 50 mm | max 8.5° | 2x110 W | 20 km ³ | | | |

1/ depending on seat adjustment

2/ total weight is depending on configuration

3/ Please note that range quoted is based upon ISO 7176-4:2008 and may vary dependent upon many factors (user weight, terrain, battery condition etc)

Please note colours may vary slightly from those displayed above.

Frame colours



Black



Transparent Purple

More detailed information, technical data, references for transport or description of certain options can be found in the wheelchair's user manual.

Invacare Limited
Pencoed Technology Park - Pencoed -
Bridgend CF35 5AQ - United Kingdom
Tel: +44 1 656 776 222
Fax: +44 1656 776 220
E-mail: uk@invacare.com
Sales Order E-mail:
ordersuk@invacare.com
www.invacare.co.uk

Invacare Ireland Ltd
Unit 5 - Seatown Business Campus
Seatown Road - SWORDS - County
Dublin - Ireland
Tel. +353 1 8107084
Fax +353 1 8107085
Email: Ireland@invacare.com
www.invacare.ie

© 2013 Invacare International Sàrl
All rights reserved.

All information quoted is believed to be correct at time of print.
Invacare reserves the right to alter product specifications without prior consultation.

Phoenix Plus - UK - 04/2013
Q99642



Yes, you can.®