






Hula

The flexible powerchair that revolves around you



Excellent indoor manoeuvrability

- The Quickie Hula's true mid-wheel drive position means it can literally turn on the spot for great indoor manoeuvrability.
- A narrow width of just 63cm (25") enables easy access through doors and confined spaces.
- A very low seat height provides easy transfers and access under tables. Seat height is adjustable between 45cm - 47cm (17 - 19") via the seat post and seat angle.

Impressive outdoor performance

- The advantages of the Quickie Hula's true mid wheel drive can be seen outdoors too, the central drive wheel increases traction and stability, and reduces sliding.
- The Quickie Hula's unique suspension system means it can climb 5cm (2") thresholds and tackle slopes of up to 6° - ideal for use in the garden or for short trips out of the home.
- Crash tested and approved to ISO 7176-19 to give you peace of mind that you're using a safe chair.



Flexible seating - now and in the future

The Hula uses the same seat as the Quickie Salsa and can be easily adjusted to a size to suit you. Choose from a wide range of options to enhance your comfort and meet your clinical needs. The Hula can be tailored to meet your requirements now and can be easily adapted for the future.



- Seat width 41 to 51cm (16 - 20")
- Seat depth 41 to 51cm (16 - 20")
- Seat height 41.5 to 51cm (16 - 20")
- Seat angle 0°, 3°, 9°
- Backrest angle 0° -12°




Impressive standard specification

- Tension adjustable backrest
- Central flip-up footplate
- Flip back armrests
- 25cm (10") drive wheels
- 20cm (8") castors
- 36Ah batteries
- Suspension
- VR2 programmable controller



Tension adjustable backrest



Flip back armrests



Central flip-up footrest

The powerchair that revolves around you

As well as the impressive standard specification the Hula is available with a variety of extras to enhance your comfort. Select from the range of JAY Comfort Seating, or more clinical options, like lateral supports - the Hula has it all.

All options are modular which means they can be easily added or removed at a later date, adding to the flexibility that's been designed into the Hula.



Manual elevating legrests



Adjustable lateral supports:
for sling or Comfort back



Available with full range of
JAY Comfort backrests
and cushions



Hula

The flexible powerchair that revolves around you



Technical Data: Seating

Seat width:	41, 43, 46, 48, 51cm (Adjustable from 41 to 51cm)	16, 17, 18, 19, 20"
Seat depth:	41, 43, 46, 48, 51cm (Adjustable from 41 to 51cm)	16, 17, 18, 19, 20"
Seat angle:	Adj. 0°, 3°, 6°	
Seat height:	45cm / 47cm (with 3° seat angle)	17" / 19" (with 3° seat angle)
Backrest angle:	Adj. -4°, 0°, 4°, 8°	
Backrest height:	Adj. 46, 48 and 51cm (sling) 57cm Comfort backrest	18, 19 and 20" 23" Comfort backrest
Backrest width:	Adjustable from 41 to 51cm	16 - 20"
Lower leg length:	34cm - 50cm	13 - 20"
Legrest angle:	90° Fixed centre mount, 80° swing-away, ELR	
Armrest height:	Adj. from 22 to 31cm	8.7 to 12"
Armrest depth:	Adj. from 31 to 41cm	12 to 16"

Specification

Overall width:	63cm	25"
Overall length: (including legrests)	97cm - 105cm	38 - 41"
Chair weight*:	82kg	13st
Turning circle:	115cm	45"
Battery capacity	36 Ah	
Range:	18km	11miles
Speed:	3 or 6kph	2 or 4 miles
Obstacle climbing ability:	5cm	2"
Maximum user weight:	136kg	21.5 stone
Max safe slope:	6°	
Ground clearance:	5cm	2"
Transit tested to ISO 7176 Part 19 and passed:	Yes	
Warranty:	2 years (please refer to owners manual for terms and conditions)	

Find out more today!

Visit www.SunriseMedical.co.uk/Hula for more information, to find your local dealer or to request a demonstration.



Get SOCIAL with QUICKIE

Keep up-to-date with the latest news on the Quickie Hula and other Quickie products through our Facebook page:

www.facebook.com/QuickieWheelchairsUK



Sunrise Medical Limited
Thorns Road
Brierley Hill
West Midlands
DY5 2LD

Tel: +44 (0) 845 605 66 88
Fax: +44 (0) 845 605 66 89

www.SunriseMedical.co.uk